

# Quaternary Ammonium Sanitizer Testing Procedure

- Know target Quat level.
- Choose correct Quat test kit.

Dip paper in quat solution, NOT FOAM SURFACE, for 10 seconds. Don't shake. Compare colors at once. Test solution should be between 65-75°F. Test solution should have a neutral pH.



**Tip:** Rinse bottles and buckets before refilling to insure proper dilution. Residual Quaternary Ammonium can elevate sanitizer level.

MICRO ESSENTIAL



[www.MicroEssentialLab.com](http://www.MicroEssentialLab.com)